

NAME: ..... ADM NO: .....

SCHOOL : .....DATE : .....

CANDIDATE'S SIGNATURE.....

233/3

CHEMISTRY

PAPER 3

TIME: 2<sup>1</sup>/<sub>4</sub>HOURS

**Kenya Certificate of Secondary Education (K.C.S.E)**

**INSTRUCTIONS TO CANDIDATES:**

- (i) Write your **name** and **index number** in the spaces provided **above**.
- (ii) **Sign** and write the **date** of examination in the spaces provided **above**.
- (iii) Answer **ALL** the questions in the spaces provided in the question paper.
- (iv) You are not allowed to start working with apparatus for the first 15 minutes of 2<sup>1</sup>/<sub>4</sub>Hours allowed for this paper. This time is to enable you read the question and make sure you have all the chemicals and apparatus required.
- (iv) Mathematical tables and silent electronic calculators **may be** used.
- (v) All working **must be** clearly shown where necessary.
- (vi) Candidates should check the question paper to ascertain that all the pages are printed as indicated and that no questions are missing

**FOR EXAMINER'S USE ONLY:**

Question	Maximum Score	Candidate's Score
1	20	
2	13	
3	7	
Total Score	40	

### **QUESTION 1**

(a) You are provided with the following

(i) Sulphuric(VI) acid labeled a solution D

(ii) Solution R prepared by dissolving 40g of NaOH in 400cm<sup>3</sup> of distilled water and made up to 1 litre.

You are required to determine the molarity of Sulphuric (VI) acid.

### **PROCEDURE**

(i) Fill the burette with Sulphuric (VI) acid.

(ii) Using pipette and pipette filter place 25cm<sup>3</sup> of solution R into conical flask.

(iii) Add 2-3 drops of phenolphthalein indicator.

(iv) Titrate solution R against solution D.

(v) Repeat the titrations to obtain two concordant titrates

(vi) Record your results in the table below.

	I	II	III
Final burette reading (cm <sup>3</sup> )			
Initial burette reading (cm <sup>3</sup> )			
Volume of solution D used (cm <sup>3</sup> )			

i) Calculate the average volume of acid solution D used . (4mks)  
(1mk)

ii) Calculate the number of moles NaOH solution (Solution R) used. (1mk)

iii) Calculate the number of moles of acid used. (1mk)

iv) Determine the molarity of Sulphuric (vi) acid. (1mk)

b) You are provided with the following;

- (i) 0.85M HCL labeled solution N.
- (ii) Sodium hydroxide labeled solution K.

You are required to determine the molar heat of neutralization **of** solution N

Procedure


- i. Measure 50cm<sup>3</sup> of solution N and transfer it into a 250cm plastic beaker provided.
- ii. Record the initial temperature of solution N to the nearest 0.5<sup>0</sup>c.
- iii. Rinse the burette thoroughly with distilled water and fill it up to the 0.0cm<sup>3</sup> mark with sodium hydroxide solution ( Solution K)
- iv. From the burette add 10cm<sup>3</sup> of solution K to solution N in the beaker. Stir gently with the thermometer and record the new temperature in the table below.
- v. Continue adding 10cm<sup>3</sup> portions of Solution K recording new temperature after each addition until 80cm<sup>3</sup> of K has been added

Volume of Solution K added (cm <sup>3</sup> )	0	10	20	30	40	50	60	70	80
Temperature <sup>0</sup> C									

(3mks)

- a) Plot a graph of temperature (y-axis) against total volume of sodium hydroxide (Solution K) (x-axis)

(3mks)



- i) From the graph determine the maximum temperature reached. (1mk)
- ii) What is the temperature change? (1mk)
- b) From the graph determine the volume of sodium hydroxide required for complete neutralization of hydrochloric acid. (1mk)
- c) Write an Ionic equation for the neutralization reaction. (1mk)
- d) Calculate the heat evolved when volume of sodium hydroxide in (c) was neutralized.  
Density of solution  $1\text{g/cm}^3$ ,  $C=4.2\text{J/g/K}$ . (1mk)
- e) Calculate the molar heat of neutralization. (1mk)

## **QUESTION 2**

You are provided with solid D. Carry out the following tests and record your observation and inferences in the spaces provided.

- a) Describe the appearance of Solid D. (1mk)

- b) Place about half of solid D in dry test. Heat it strongly and test for any gas produced using a blue and red litmus paper.

Observation	Inferences
1mk	1mk

- c) Place the rest of Solid D in boiling tube and add about 10cm<sup>3</sup> of distilled water. Shake well and add about 2cm<sup>3</sup> portions for each of the test below.

Observation	Inferences
1mk	1mk

- i) To one portion, add aqueous NaOH drop wise until in excess

Observation	Inferences
1mk	$\frac{1}{2}$ mk

ii) To the second portion, add aqueous ammonia drop wise until in excess

Observation	Inferences
1mk	$\frac{1}{2}$ mk

iii) To a third portion add about 5 drops of sodium chloride solution.

Observation	Inferences
1mk	$\frac{1}{2}$ mk

iv) To a fourth portion, add dilute Barium nitrate solution.

Observation	Inferences
1mk	$\frac{1}{2}$ mk



Put the remaining portion of solid J in a boiling tube and add about 8cm<sup>3</sup> of distilled water. Shake to dissolve.

- ii) To about 2cm<sup>3</sup> of solution J in a test tube add 2 to 3 drops of bromine water.

Observation	Inferences
1mk	1mk

- iii) To about 2cm<sup>3</sup> of solution J in a test tube, add about 1cm<sup>3</sup> of acidified potassium dichromate (vi). Warm gently and allow to stand for a minute.

Observation	Inferences
1mk	$\frac{1}{2}$ mk

- iv) To about 2cm<sup>3</sup> of solution in a test tube, add a small amount of solid Sodium hydrogen carbonate.

Observation	Inferences



1mk	$\frac{1}{2}$ mk
-----	------------------