





# MARKING SCHEME, PRACTICAL TERM 1 2024

233/3

## CHEMISTRY PAPER

3 (Practical)

2¼ Hours

### INSTRUCTIONS:

- (i) Write your name and index Number in the spaces provided.
- (ii) Answer *All* questions in the spaces provided in the question paper.
- (iii) You are *NOT* allowed to start working with the apparatus for the first *15 minutes of the 2 ¼ hours* allowed for this paper. This time is to enable you to read the question paper and make sure you have all the chemicals and apparatus that you may need.
- (iv) All working *MUST* be clearly shown where necessary.
- (v) Mathematical tables and electronic calculators may be used.

### FORM EXAMINER'S USE ONLY

Question	Max. Score	Score
1	20	
2	14	
3	06	
Total Score	40	

## QUESTION 1

You are provided with:

- Solid Z ; 5.0 g of  $(\text{COOH})_2 \cdot X \text{H}_2\text{O}$
- Solution Y ; 0.125658 M  $\text{KMnO}_4$

You are required to :

- a) Determine the solubility of Z at different temperatures.
- b) Determine the number of moles of water of crystallization in solid Z.

## PROCEDURE 1

- a) Using a burette, add  $4\text{cm}^3$  of distilled water to solid Z in a boiling tube.
  - Heat the mixture while stirring with the thermometer to about  $80^\circ\text{C}$ .
  - When the whole solid dissolves, allow the solution to cool while stirring with the thermometer
  - Note the temperature at which crystals first appear and record this temperature in the table 1 below.
- b) Using a burette add  $2\text{cm}^3$  more distilled water into the content of the boiling tube and warm until the solid dissolve.
  - Remove from the flame and allow the solution to cool in air while stirring.
  - Record the temperature at which crystal first appears in table 1.
  - Repeat procedure (b) 3 more times and complete table 1 below.
  - Retain the content of the boiling tube for procedure II

Table I.....,..... 6mks

Volume of water in the boiling tube ( $\text{cm}^3$ )	Temperature at which crystals of solid A appear ( $^\circ\text{C}$ )	Solubility of solid Z (g/100g of water)
-------------------------------------------------------	----------------------------------------------------------------------	-----------------------------------------

4		
6		
8		
10		
12		

I). COMPLETE TABLE..... 3MKS

- all temp 2mk, 4 temp readings 1 ½ mk, 3 temp readings 1mk, 2 temp ½ mk
- all solubility well calculated 1mk, 3-4 solubility calculated well ½ mk
- penalize ½ mk ONCE for temp reading above 80 and below 30
- award a max 1 ½ mk for constant temp readings

II. DECIMALS... 1mk


- TIED to all temperatures either constant or 1 dp of .0 or .5  
Otherwise award 0 mk

III. ACCURACY..... 1mk

Tied to the 1<sup>st</sup> temp reading , if within + or – 2<sup>0</sup> of the school value ,award 1 mk, otherwise award 0 mk

IV. TREND..... 1mk continuous drop in temperature, otherwise award 0 mk

I. a) Draw a graph of solubility of solid Z (vertical axis) against temperature  
(3mks)



- LABELLING.....  $\frac{1}{2}$  mk  
both axes labelled for  $\frac{1}{2}$  mk. penalize fully for inverted axes. penalize fully for wrong units otherwise ignore if units are not given
- SCALE.....  $\frac{1}{2}$  mk ( check the space for the actual plots and if  $\frac{1}{2}$  of the boxes are used in both axes, award fully), the scale should accomodate all plots, the spacing must be consistent.
- PLOTTING.... 1 mk ( 5 plots correctly plotted.....1MK, 3-4 correctly plotted  $\frac{1}{2}$  mk)
- CURVE..... 1mk ( a smooth curve passing through at least 3 plots which include the 1<sup>st</sup> plot , the last plot and any other one plot). The curve should not be extrapolated

b) From your graph determine the solubility of solid Z at 60°C  
(1mk)

- correct reading from your graph

## PROCEDURE II

a) – Transfer the contents of the boiling tube into a 250ml volumetric flask.

- Add distilled water up to the mark

- Label this solution Z

b) – Using a clean pipette and a pipette filler, transfer 25ml of solution Z into a conical flask.



- Warm the mixture up to  $60^{\circ}\text{C}$
  - Fill a burette with solution Y
  - Titrate Y against the hot solution Z until a permanent pink colour persist
  - Record your results in Table 2 below
- c) Repeat (b) 2 more times and record your results in the table 2 below.

TABLE 2

	I	II	III
FINAL BURETTE READING			
INITIAL BURETTE READING			
VOLUME OF SOLUTION Y USED ( $\text{cm}^3$ )			

(4mk)

TABLE 2

Complete table -----1mk

- (i) Complete the table with 3 titrations done --(1mk)
- (ii) Incomplete table with 2 titrations done --- (1/2mk)
- (iii) Incomplete table with 1 titration done --(0mk)

#### Penalties

- I. Wrong arithmetic
- II. Inverted table
- III. Unrealistic value

Penalize  $\frac{1}{2}$  mark for each to a maximum of  $\frac{1}{2}$ mk

Decimals -----1mk

(Tied to the first and second row only )

Conditions

Accept either 1 or 2 decimal points used consistently

If the 2<sup>nd</sup> decimal point is used . can only be 0 or 5



Accuracy -----1mk

Compare any correct titre value in the 3<sup>rd</sup> row with the school value (sv)

Conditions

I.If within  $\pm 0.1\text{cm}^3$  of sv ---1mk

II.If not within  $\pm 0.1\text{cm}^3$  but within  $\pm 0.2\text{cm}^3$  of sv ---1/2mk

III.Beyond  $\pm 0.2\text{cm}^3$  of sv -----0mk

NB/ if there is wrong arithmetic in the table compare the sv with the correct value and credit accordingly.

Principle of averaging ----1mk

Value average must be shown and must be within  $\pm 0.2\text{cm}^3$  of each other; conditions.

I.3 values averaged and consistent -1mk

II.3 values done and only 2 possible averaged -1mk

III.2 titrations done and averaged -1mk

IV.2 titrations done and inconsistent -0mk

V.3 titrations done and consistent but only two averaged -0mr

Final accuracy -1mk

Compare correct student average titre with the (sv )

I.If within  $\pm 0.1$  of sv -1mk

II.If within  $\pm 0.2$  of sv -1/2mk

III.If beyond  $\pm 0.2$  of sv -0mk

NB// If the candidate has averaged wrong values pick the correct value if any ,average and credit accordingly

II) a) Calculate the average volume of solution Y used (1mk)

Marked in the table

b) Calculate the number of moles of Y used 1mk

$$\frac{\text{aver volume}}{1000} \times 0.125658 \frac{1}{2} \text{ mk} = \text{corr ans } \frac{1}{2} \text{ mk}$$

c) Given 2 moles of  $\text{KMnO}_4$  react with 5 moles of Z, calculate the number of moles of Z in  $25\text{cm}^3$  (1mk)

mole ratio 2 : 5

$$\text{moles of A} = \frac{\text{answer in b aboven}}{2} \times 5 \frac{1}{2} \text{ mk} = \text{corr ans } \frac{1}{2} \text{ mk}$$



d) Calculate the molarity of Z (1mk)

$$= \frac{\text{answ in c}}{25} \times 1000 \text{ } \frac{1}{2} \text{ mk} = \text{corr ans } \frac{1}{2} \text{ mk}$$

e) Determine the molar mass of Z (1mk)

$$= \frac{5}{\text{ans in d}} \text{ } \frac{1}{2} \text{ mk} = \text{corr ans } \frac{1}{2} \text{ mk}$$

f) Determine the value of X (C=12, O=16 H=1) (1mk)

$$90 + 18x = \text{ans in e } \frac{1}{2} \text{ mk}$$

$$X = \frac{\text{ans in e} - 90}{18} = \text{corr ans } \frac{1}{2} \text{ mk}$$

2. You are provided with solid W. Carry out the tests below and record your observations and inferences in the spaces provided.

a) Place solid W in a dry boiling tube. Add about 20cm<sup>3</sup> distilled water and shake well. Filter to obtain the residue and the filtrate.

Observations (1mk)	Inferences (1mk)
<ul style="list-style-type: none"> <li>- colourless filtrate <math>\frac{1}{2}</math> mk</li> <li>- white residue <math>\frac{1}{2}</math> mk</li> </ul>	<ul style="list-style-type: none"> <li>- is a mixture of a soluble and insoluble salt</li> </ul>

b) Divide the filtrate into 3 portions.

i) To the first portion, add Ba (NO<sub>3</sub>)<sub>2</sub>, followed by 5 drops of dilute nitric (v) acid

Observations ( 1mk)	Inferences ( $\frac{1}{2}$ mk)
<ul style="list-style-type: none"> <li>- white ppt <math>\frac{1}{2}</math> mk insoluble in the acid <math>\frac{1}{2}</math> mk</li> </ul>	SO <sub>4</sub> <sup>2-</sup>

ii) To the 2<sup>nd</sup> portion, add aqueous ammonia drop wise until in excess

Observations ( 1mk)	Inferences (1mk)
<ul style="list-style-type: none"> <li>- white ppt <math>\frac{1}{2}</math> mk insoluble in excess <math>\frac{1}{2}</math> mk</li> </ul>	-Mg <sup>2+</sup> 1/2 mk, Al <sup>3+</sup> present 1/2 mk Penalize $\frac{1}{2}$ mk for @ contradictory ion including Pb <sup>2+</sup>

iii) To the 3<sup>rd</sup> portion, add aqueous sodium hydroxide dropwise until in excess



Observations ( 1mk)	Inferences ( ½ mk)
- white ppt ½ mk insoluble ½ mk in excess	- $Mg^{2+}$

c) Wash the residue and put it in a test tube. Add about  $15\text{cm}^3$  of dilute nitric (v) acid. Test the gas produced with a burning splint. Divide the solution into 3 portions

Observations ( 1mk)	Inferences (1mk)
<ul style="list-style-type: none"> <li>- bubbles ½ mk of a colourless gas</li> <li>- put off a burning splint ½ mk</li> <li>- residue dissolves ½ mk</li> </ul>	$\text{CO}_3^{2-}$ present 1mk

i)Add sodium hydroxide to the first portion dropwise until in excess to the 1<sup>st</sup> portion

Observations ( 1mk)	Inferences (1mk)
White ppt ½ mk soluble ½ mk	$\text{Zn}^{2+}$ , $\text{Pb}^{2+}$ , $\text{Al}^{3+}$ present

ii)Add 2 drops of potassium iodide to the 2<sup>nd</sup> portion.

Observations ( ½ mk)	Inferences ( ½ mk)
Yellow ppt ½ mk	$\text{Pb}^{2+}$

iii)Add  $2\text{cm}^3$  of solution R to the 3<sup>rd</sup> portion and wait for few seconds.

Observations ( 1mk)	Inferences (1mk)
White ppt 1mk	Soln R contains $\text{Cl}^-$ , $\text{SO}_4^{2-}$ , $\text{SO}_3^{2-}$ , $\text{CO}_3^{2-}$ , or $\text{Br}^-$ ions Any 4 ions 1mk Any 3 ions ½ mk



3. You are provided with solid V. Carry out the tests below and record your observations and inferences in the spaces provided.

a) Place a half spatula of solid V in a dry boiling tube. Add about 6cm<sup>3</sup> distilled water and shake well.

Observations (1mk)	Inferences ( ½ mk)
- dissolves ½ to form a colourless solution ½	Polar organic compound 1mk

b) Divide into 3 portions

i) To the first portion, add ½ spatula of NaHCO<sub>3</sub>

Observations ( ½ mk)	Inferences (1/2 mk)
No bubbles	RCOOH or H <sup>+</sup> absent

ii) To the 2<sup>nd</sup> portion, add 3 drops of acidified K<sub>2</sub>Cr<sub>2</sub>O<sub>7</sub> and warm

Observations ( ½ mk)	Inferences ( ½ mk)
- acidified K <sub>2</sub> Cr <sub>2</sub> O <sub>7</sub> remain orange	ROH absent

iii) To the 3<sup>rd</sup> portion, add 3 drops of acidified KMnO<sub>4</sub>

Observations ( 1mk)	Inferences ( ½ mk)
Acidified KMnO <sub>4</sub> remain purple	=C=C= absent

b) Scoop the remaining with a spatula and heat it on a non luminous flame

Observations ( ½ mk)	Inferences ( ½ mk)
-yellow sooty flame	Long chain saturated organic cpd ½ mk -ignore mention of unsaturated organic compound present or =C=C= present

